## Hydrophile Acryl Intraokularlinsen i-Flex-6025C-6011C-i-Soft-7011C-6513C

### PRODUKTBESCHREIBUNG

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## Hydrophilic acrylic intraocular lenses i-Flex-6025C-6011C-i-Soft-7011C-6513C

### PRODUCT DESCRIPTION

**i-Medical**<sup>®</sup> intraocular lenses are mono-focal folding acrylic IOLs. These intraocular lenses also have a UV-blocker to reduce the high-energy light in the ultraviolet light spectrum.They are implanted in the lens capsule behind the iris in order to take over the vision function of the natural crystalline lens.

**i-Medical**<sup>®</sup> intraocular lenses are manufactured to the sharp-edge design for the reduction or prevention of clouding of the rear capsule with cell proliferations.

**i-Medical**<sup>®</sup> intraocular lenses (IOLs) of hydrophilic acrylic are optical implants made of copolymers of the highest medical quality. The water content is 26 %. The lens has a high refractive index of 1.5015 in the dry condition and 1.460 in the moist condition. Hydrophilic acrylic is a high-grade bio-compatible material.

The hydrophilic acrylic copolymer material has a refractive index of 1.46 and a transmission of visible light of > 90 %. The cut-off wavelength of 10 % transmission for UV light is 376 nm for a +4 diopter lens and 390 nm for a +34 diopter lens.

**i-Medical**<sup>®</sup> IOLs are made from a single piece, ground and surface-finis hed, and are intended as surgical implants to replace the human lens. The extraction can be carried out extra-capsular and / or by phakoemulsification. The IOL of hydrophilic acrylic is available from 0.0 to 30.0 dioptres in steps of 0.5. The optical quality and the dioptre strength of every single IOL are checked individually.

**INDICATIONS APHAKIA** after cataract surgery

**DETERMINATION OF THE DIOPTR POWER TO BE IMPLANTED**
It can be determined by methods described in published literature.

### CONTRAINDICATIONS

Apart from non-specific contraindications related to any form of ocular surgery, the following specific contraindications must be respected (Non-exhaustive list):

- Bilateral congenital cataract
- Microphthalmia
- Chronic uveitis
- Active ocular diseases (active diabetic retinopathy), glaucoma
- Patients receiving or having received chloroquine (Nivaquine<sup>™</sup>) treatment.
- Children under the age of 12 months
- Corneal dystrophy or endothelial insufficiency

### COMPLICATIONS RELATED TO IOL IMPLANTATION

The possible complications related to IOL implantation have been listed according to those described in literature (Non-exhaustive list):

Pre-operative complications:

- Anterior vitrectomy
- Excessive vitreous loss
- Expulsive hemorrhage
- Posterior capsule rupture
- Iris damage

### Post-operative complications:

- Elevated IOP
- Cystoid macular edema opacification
- Retinal detachment
- Atrophy of the Iris
- Anterior segment inflammation
- Synechia
- Anterior chamber collapse
- IOL Precipitates
- Endophthalmitis
- IOL opacification
- Posterior capsule opacification
- Hypopyon
- High anisotropia and aniseikonia
- Corneal edema
- Pupillary block
- IOL-decentering
- IOL dislocation and subluxation

#### STERILIZATION AND PACKAGING

This implant is sterilized by humid heat before being placed in an indivi-dual sterile blister. A patient card and self-adhesive labels identifying the implant facilitate for the patient’s medical follow-up are included

### INSTRUCTIONS

The cardboard box is opened by pulling on the flap indicated. Remove the sterile package and peel open smoothly and take the IOL blister out. Open the blister smoothly and grasp the lens using non-toothed forceps. All manipulations required for implantation must be performed using non-toothed forceps and polished instruments.

### PRECAUTIONS:

**DO NOT RESTERILIZE – SINGLE USE ONLY!**

After opening the blister bag, the IOL must be used immediately. Check the integrity of the external sterile packaging before use.

### WARNING

**DO NOT USE IF THE EXTERNAL STERILE PACKAGING IS DAMAGED!**
The surgeon is informed that, because of its hydrophilicity, the IOL can absorb substance with which it is placed in contact (e.g. antibiotics, viscoelastic).
**TRANSFERTING FROM COLD TEMPERATURES TO WARM TEMPERA-TURES MAY CAUSE THE LENS TO CLOUD. DO NOT OPEN OR USE UNTIL THE LENS REACHES THE ROOM TEMPERATURE AND BECOMES CLEAR. DO NOT USE AFTER THE EXPIRY DATE.**
When reused, there is an increased RISK for material damage and increased risk of infection.
**For Single Use Only!**

### STORAGE CONDITIONS

Store in a cool dry place (max 45° C, min 0° C).

**i-Medical**<sup>®</sup> declines all responsibility following the use of an implant or related to the development of pre and post-operative complications.

## Lentes intraoculares de acrílico hidrófilo i-Flex-6025C-6011C-i-Soft-7011C-6513C

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