Faltbare Hydrophobe Intraokularlinse

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Foldable Hydrophobic Intraocular Lens

OVERVIEW AS natural – hydrophobic, aspheric – Preloaded IOL

INDICATIONS

The **i-Medical®** OVERVIEW AS natural (aspheric) hydrophobic foldable intraocular lenses are optical implants for the replacement of the human lens for visual correction of aphakia / cataracts. The lenses are intended to be placed in the capsular bag only.

DESCRIPTION

i-Medical® OVERVIEW AS natural hydrophobic one-piece IOL, manufactured from a hydrophobic acrylate co-polymer.

Optic Material: hydrophobic acrylate co-polymer, with natural blue light filter.
 IOL Design: Modified C Loop
 Biconvex Optic, aspheric
 Diopter range: +10.0 to +30 diopter in 0.5 diopter steps
 Sterilisation type: EO sterilisation
 The hydrophobic acrylic copolymer material has a refractive index of 1.525 and a transmission of visible light of > 90 %. The cut-off wavelength of 10 % transmission for UV light is 376 nm for a +4 diopter lens and 390 nm for a +34 diopter lens.

The **i-Medical®** Overview AS natural hydrophobic technology:

- The **i-Medical®** OVERVIEW AS natural hydrophobic one-piece IOLs are manufactured using the cast molding process with 360 degree square edge technology.
- **i-Medical®** OVERVIEW AS natural hydrophobic one-piece IOLs have a special manufacturing technology which provides a natural blue light filter – like the human eye at the age of about 40 years.

- **i-Medical®** OVERVIEW AS natural is able to filter high energy UV light harmful to the human eye.
- **i-Medical®** OVERVIEW AS natural hydrophobic one-piece IOLs have an aberration-free neutral aspheric optic.
- The **i-Medical®** OVERVIEW AS natural hydrophobic one-piece IOLs are "glistering-free" IOLs.
- **i-Medical®** OVERVIEW AS natural hydrophobic one-piece IOLs have been developed in order to minimise PCO.

Implantation recommendation of the i-Medical® Preloaded IOL Overview AS natural:

See the enclosed **i-Medical®** recommendation / instruction Steps 1 to 9.

The following indications and contraindications are based on research of medical literature and are to be used only as guides. The list is indicative and makes no claim to completeness.

CONTRAINDICATIONS

Implantation is not advisable when the IOL may aggravate an existing condition, interfere with the diagnosis or the treatment of pathology or present a risk to the sight of the patient. Such conditions include:

1. Chronic or recurrent uveitis
2. Proliferative diabetic retinopathy
3. Corneal endothelial dystrophy
4. Acute eye disease or infection, external or internal
5. Severe complication during surgery
6. Choroidal vitreous loss
7. Choroidal haemorrhage
8. Non-age-related cataract
9. Microphthalmus
10. Suspected microbial infection
11. Medically uncontrolled glaucoma
12. Severe optic atrophy
13. Uncontrollable glaucoma

Complications related with intraocular lens implantation:

1. Capsular rupture
2. Vitreous loss
3. Lens decentration and luxation

WARNINGS

Among those directly linked to the IOL are: Decentration, subluxation, precipitates on the IOL surface. As with all surgical procedures, cataract surgery with IOL implantation presents risks, which the surgeon must consider and evaluate the risk/benefit ratio. Some of the potential complications of cataract surgery are:

1. Wound leak
2. Retinal detachment
3. Cystoid macular oedema
4. Corneal decompensation
5. Corneal oedema
6. Pupillary block
7. Iritis
8. Endophthalmitis
9. Secondary glaucoma

Some of the listed complications may require a second surgical intervention. The complications may occur following implantation of any intraocular lens.

PRECAUTIONS

1. An IOL should only be implanted in the human eye by an experienced ophthalmic surgeon following corresponding indication.
2. Before the opening of the blister packaging, the IOL configuration, power and expiration date should be verified before opening the blister pack.
3. Do not use the intraocular lens if the sterile package has been opened or damaged or if the expiration date has already been exceeded.
4. Do not restylize the lens.
5. The intraocular lens must be implanted only in the capsular bag of the human eye.

DIRECTIONS FOR USE

Open the sterile **i-Medical®** blister, and remove the blister closure. Take out the sterile injector with the preloaded **i-Medical®** intraocular lens from the blister. Pull and remove the yellow protective cover in the direction of the arrow. Apply 5-6 drops of BSS solution and also viscoelastic fluid i-Visc[®] onto the lens surface of the **i-Medical®** IOL in the cartridge. Close the flap of the cartridge until it engages. Now fill the injection opening of the cartridge with viscoelastic fluid i-Visc[®]. Now push the **i-Medical®** IOL with the piston of the injector into the 1st quarter of the of the injection opening of the cartridge. Ensure that it slides smoothly and evenly. The preloaded **i-Medical®** intraocular lens is now "ready to use" and can be implanted into the eye.

When used, there is an increased RISK for material damage and increased risk of infection. For Single Use Only!

STORAGE

DO not expose the lens to direct sunlight or temperatures above 45 °C. Do not store at temperatures below 0 °C.

Implant the intraocular lens only at the normal room temperature in the operating room.

HOW SUPPLIED

i-Medical® OVERVIEW AS natural is supplied sterile in the injection system packed in a blister pack. The lenses are carefully checked and inspected by the manufacturer to ensure a high product quality. If a defect or deformation is noted or suspected, the lens should be returned to the manufacturer.

A-KONSTANT

The „A“ constant value recommended by the manufacturer is specified on the outer packaging of the IOL. This number is given only as a guideline and is not based on exact clinical data. It is advised that the surgeon calculates his own value for each type and specification of intraocular lens.

CALCULATION OF LENS POWER

For every patient to be implanted with an intraocular lens, the ophthalmic surgeon must calculate the exact refractive power (dpt) of the lens to be implanted and select the IOL configuration suitable for the patient accordingly.

LIABILITY DISCLAIMER

- i-Medical®** will not be liable for any injury suffered by the patient as a result of:
 - Any implantation method or technique used by a physician to implant the lens which is not recommended here by **i-Medical®**.
 - The selection of the IOL configuration and refractive power on the intraocular lens by the ophthalmic surgeon.
 - For operating hygiene in the area of the OP of the ophthalmic surgeon implanting the IOL.

i-Medical® guarantees only that the IOL product in question complies with the product properties assured here and the standards of the relevant European CE approval for medical products. The success of the implantation of this IOL for the patient is the sole responsibility of the ophthalmic surgeon performing the operation, for which **i-Medical®** accepts no liability.

REPORTING

Adverse reactions and potentially sight-threatening complications that may reasonably be regarded as related to the intraocular lens must be reported immediately and in writing to **i-Medical®**. This information is requested from all the implanting surgeons in order to document potentially long-term effects of IOL implants. We can thereby continue to ensure the patient and product safety in the long term.

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