



BioLine® YELLOW



ACCURATE ASPHERIC

Type	ACCURATE® ASPHERIC One-piece hydrophilic acrylic IOL UV / Blue light filter
Optic size	6.0 mm
Overall length	12.0 mm
Haptic angulation	0°
Optic	Biconvex / Aspherical
Central thickness	0.97 mm (+20.0 Dpt.)
Design	360° square edges
Material	Hydrophilic acrylic copolymer 26% water
Available diopters	-20.0 Dpt. to +60.0 Dpt. in 0.5 Dpt. increments „Diopter below -10.00D and diopters higher than 40D are manufactured only on special request.“
Sterilization	Steam autoclaved
Refractive index	1.46
Estimated A-Constant	118.8
Estimated anterior chamber depth	4.98 mm
Recommended incision size	2.2 / 2.5 mm
SA	0.0 neutral
ABBE	58
Injector/Cartridge	Attached

INTRAOCCULAR LENSES SQUARE EDGE TECHNOLOGY

