



BioLine® YELLOW

# BLUELIGHT

<b>Type</b>	BLUELIGHT One-Piece hydrophilic acrylic IOL UV/Blue light filter
<b>Optic size</b>	6.0 mm
<b>Overall length</b>	12.5 mm
<b>Haptic design</b>	Modified C
<b>Haptic angulation</b>	0°
<b>Optic</b>	Biconvex / Spherical
<b>Central thickness</b>	0.98 mm (+20.0 Dpt.)
<b>Design</b>	360° square edges
<b>Material</b>	Hydrophilic acrylic copolymer 26% water
<b>Available diopters</b>	-20.0 Dpt. to +40.0 Dpt. in 0.5 Dpt. increments, Diopters below -1.0 and above +30.0 are manu- factured only on special request
<b>ABBE</b>	58
<b>Sterilization</b>	Steam autoclaved
<b>Refractive index</b>	1.46
<b>Estimated A-Constant</b>	118.0
<b>Anterior chamber depth</b>	4.96 mm
<b>Recommend incision size</b>	2.2 / 2.5 mm
<b>SA neutral</b>	0.0
<b>i-Ject Set (Injector + Cartridge)</b>	Size small or medium

## INTRAOCCULAR LENSES SQUARE EDGE TECHNOLOGY

